FPA anticoagulant solution

REF ARO13A

ACS Vial of 20 mL / Flacon de 20 mL

FPA anticoagulant solution
FOR RESEARCH USE ONLY.
NOT FOR USE IN DIAGNOSTIC PROCEDURES

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English, revision: 12-2023

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INTENDED USE:

Anticoagulant solution for use in Fibrino-Peptide A (FPA) assay.
This kit is for research use only and should not be used for patient diagnosis or treatment.

REAGENTS:

ACS FPA anticoagulant solution : special anticoagulant solution for FPA assay, solution containing TriSodium citrate, heparin, aprotinin, hirudin, and sodium azide (0.9 g/L), liquid form. Contains preservatives and stabilizers.

The product is classified as non-hazardous and is not subject to labeling according to EC Regulation No. 1272/2008 [CLP].

WARNING AND PRECAUTIONS:

- Waste should be disposed of in accordance with applicable local regulations.
- Use only the reagents from the same batch of kits.
- This device of *in vitro* use is intended for professional use in the laboratory.

REAGENT PREPARATION:

Reagent is ready to use.

REAGENTS AND MATERIALS REQUIRED BUT NOT PROVIDED:

HYPHEN BioMed ZYMUTEST™ FPA (RK016B) reagent.

SPECIMEN COLLECTION AND PREPARATION:

The blood (9 volumes) should be carefully collected onto the **ACS** anticoagulant (1 volume) (trisodium citrate and inhibitors) by clean venipuncture. Discard the first tube. Promptly mix with **ACS** and centrifuge (eg 20 min at 2500g.) to collect the plasma supernatant thus ready for bentonite adsorption. Bentonite treatment should be promptly done, if not possible plasma should be promptly frozen at -20°C or below. Just before use, thaw plasma at 37°C and proceed to bentonite treatment. Specimen should be promptly prepared and stored in accordance with applicable local recommendations.

STORAGE AND STABILITY:

Unopened reagents should be stored at 2-8°C in their original packaging. Under these conditions, they can be used until the expiry date printed on the kit. Reagent stability after opening, free from any contamination or evaporation, and stored closed, is of 4 weeks at 2-8°C.

LIMITATIONS:

- Any reagent presenting no limpid appearance or showing signs of contamination must be rejected.
 - Refer to the specific associated kit inserts for instructions concerning the test to be performed.
- The results obtained should be for research purposes only and not used for patient diagnosis or treatment.**

eIFU is available on the HYPHEN BioMed website.

For customer support, please contact your local provider or distributor (see www.hyphen-biomed.com).

Changes compared to the previous version.

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