Imidazole buffer

- REF** ARO21K **Bf** 3 vials x 25 mL / 3 flacons x 25 mL
- REF** ARO21L **Bf** 12 vials x 25 mL / 12 flacons x 25 mL
- REF** ARO21B **Bf** 1 vial x 15 mL / 1 flacon x 15 mL
- REF** ARO21M **Bf** 3 vials x 15 mL / 3 flacons x 15 mL
- REF** ARO21N **Bf** 12 vials x 15 mL / 12 flacons x 15 mL

Imidazole buffer

Tampon Imidazole

English, revision: 03-2022

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INTENDED USE:

Dilution buffer for hemostasis assays using an automated method.
 This device of *in vitro* diagnostic use is intended for professional use in the laboratory.

REAGENTS:

Bf Solution containing Imidazole at approximately 0.05 M, pH 7.30 ± 0.20, liquid form. Contains preservatives and stabilizers.
 H360D: May damage the unborn child.

WARNING AND PRECAUTIONS:

- Waste should be disposed of in accordance with applicable local regulations.
- Use only the reagents from the same batch of kits.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

REAGENT PREPARATION:

Bf Reagent is ready to use, homogenize and load it directly on the analyzer following Application Guide instruction.
 For another use, allow to stabilize for 30 minutes at room temperature (18-25°C), homogenize before use.

STORAGE AND STABILITY:

Unopened reagents should be stored at 2-8°C in their original packaging. Under these conditions, they can be used until the expiry date printed on the kit.

Bf Reagent stability after opening, free from any contamination or evaporation, and stored closed, is of:

- 60 days at 2-8°C.
- Stability on board of the analyzer: see the specific Application Guide.

Combination of storage are not recommended.

REAGENTS AND MATERIALS REQUIRED BUT NOT PROVIDED:

- For use with HYPHEN BioMed hemostasis reagents on CS-series, STA-R® family, ACL-TOP® family, CN-series.

LIMITATIONS:

- Any reagent presenting no limp appearance or showing signs of contamination must be rejected.
- Refer to the specific associated kit inserts for instructions concerning the test to be performed.

RESULTS:

- Lot to lot variability measured on 3 lots is: %CV ≤10%.

eIFU is available on the HYPHEN BioMed website.

Changes compared to the previous version.

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