# Special Tris-BSA buffer

REF AR026A Bf 3 vials x 25 mL / 3 flacons x 25 mL  
 REF AR026K Bf 12 vials x 25 mL / 12 flacons x 25 mL

Special Tris-BSA buffer

Tampon special Tris-BSA

English, revision: 03-2022

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### INTENDED USE:

Special dilution buffer for chromogenic assays using an automated method. This device of *in vitro* diagnostic use is intended for professional use in the laboratory.

### REAGENTS:

**Bf** Tris-BSA Buffer, pH : 7.95 ± 0.15. Ready to use. Contains BSA, human AGP, preservatives and stabilizers.

The product is classified as non-hazardous and is not subject to labeling according to EC Regulation No. 1272/2008 [CLP].

### WARNING AND PRECAUTIONS:

- Some reagents provided in these kits contain materials of human and animal origin. Whenever human plasma is required for the preparation of these reagents, approved methods are used to test the plasma for the antibodies to HIV 1, HIV 2 and HCV, and for hepatitis B surface antigen, and results are found to be negative. However, no test method can offer complete assurance that infectious agents are absent. Therefore, users of reagents of these types must exercise extreme care in full compliance with safety precautions in the manipulation of these biological materials as if they were infectious.
- Waste should be disposed of in accordance with applicable local regulations.
- Use only the reagents from the same batch of kits.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

### REAGENT PREPARATION:

**Bf** Reagent is ready to use, homogenize and load it directly on the analyzer following Application Guide instruction.

### STORAGE AND STABILITY:

Unopened reagents should be stored at 2-8°C in their original packaging. Under these conditions, they can be used until the expiry date printed on the kit.

**Bf** Reagent stability after opening, free from any contamination or evaporation, and stored closed, is of:

- 60 days at 2-8°C.
- Stability on board of the analyzer: see the specific Application Guide.

Combination of storage are not recommended.

### REAGENTS AND MATERIALS REQUIRED BUT NOT PROVIDED:

- For use with HYPHEN BioMed chromogenic reagents (BIOPHEN™ FVIII (221402 / 221406), BIOPHEN™ FVIII variant (227102) and BIOPHEN™ FIX (221801 / 221802 / 221806)) on CS-series, STA-R® family, ACL-TOP® family, CN-series.

### LIMITATIONS:

- Any reagent presenting no limp appearance or showing signs of contamination must be rejected.
- Refer to the specific associated kit inserts for instructions concerning the test to be performed.

### RESULTS:

- Lot to lot variability measured on 3 lots is: %CV ≤ 10%.

eIFU is available on the HYPHEN BioMed website.

Changes compared to the previous version.

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