Tris-NaCl-EDTA buffer, pH 7.85

Tampon Tris-NaCl-EDTA, pH 7,85

1 vial of 20 mL / 1 flacon de 20 mL
Set of 4 x 20 mL / Coffret de 4 x 20 mL
REF ARO32A / ARO32K

FOR *IN VITRO* DIAGNOSTIC USE.POUR USAGE DE DIAGNOSTIC *IN VITRO*.**LOT****EXP**

English, last revision 06-2018

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INTENDED USE:

Special dilution buffer reducing Heparin Interferences for anti-Xa chromogenic assays, especially using BIOPHEN™ DiXal kit (221030) or BIOPHEN™ Heparin LRT kits (221011/221013/221015)

REAGENTS:

[Bf] Tris-NaCl-EDTA buffer, pH 7.85. Ready to use. Contains 1% PEG and an heparin neutralizing substance.
1 vial of 20 mL (AR032A).
4 vials of 20 mL (AR032K).

Reagent contains small amounts of sodium azide (0.9 g/L), see WARNINGS AND PRECAUTIONS.

WARNING AND PRECAUTIONS:

- In contact with lead or copper pipes, sodium azide can generate explosive compounds.
- Waste should be disposed of in accordance with applicable local regulations.
- Handle the reagents with care to avoid contamination during use. If possible, avoid reagent evaporation during use by limiting the liquid-air exchange surface. Evaporation reduces the reagent's stability in the analyzer.
- Aging studies show that the reagents can be shipped at room temperature without degradation.
- To preserve reagent stability, seal the vials after use with their respective caps.
- For *in vitro* diagnostic use.

REAGENT PREPARATION AND STABILITY:

[Bf] Ready to use. Allow to stabilize for 30 minutes at room temperature (18-25°C), before use.

Homogenize the reagent prior to use, while avoiding formation of foam.

Reagent stability after opening, excluding any contamination or evaporation, and stored in the original vial, is of:

- 2 months at 2-8°C.
- 7 days at room temperature (18-25°C).

STORAGE CONDITIONS:

Unopened reagents should be stored at 2-8°C in their original packaging. Under these conditions, they can be used until the expiry date printed on the kit.

LIMITATIONS:

- Any reagent presenting an unusual appearance or showing signs of contamination must be rejected.
- Refer to the specific associated kit inserts for instructions concerning the tests to be performed.

SYMBOLS:

Symbols used and signs listed in the ISO 15223-1 standard, see Symbol definitions document.

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