

## CRYOcheck™ CorPac™

*For Research Use Only. Not for use in diagnostic procedures.*

### **Product Identification**

Catalog Number: **CCCP-15**

Lot Number: **CP023**

Expiry: **2025-11-30**

Vial Number	PT (sec)	APTT (sec)	Clauss Fibrinogen (g/L)
1	13.3	35.2	3.39
2	13.2	40.3	3.24
3	12.6	32.5	3.35
4	11.8	33.0	3.03
5	12.5	36.6	4.45
6	12.6	31.4	3.00
7	12.6	34.3	4.25
8	30.7	53.5	1.48
9	33.6	59.4	1.39
10	35.5	62.7	1.33
11	37.0	66.0	1.27
12	39.8	68.8	1.18
13	42.2	74.4	1.11
14	45.8	79.5	0.99
15	49.0	87.1	0.97
16	53.6	93.3	0.83
17	58.1	104.9	0.79
18	34.2	55.6	1.62
19	37.3	58.8	1.57
20	41.0	64.1	1.41
21	43.3	67.2	1.32
22	45.6	71.3	1.22
23	47.8	74.8	1.17
24	51.4	80.0	1.11
25	55.1	86.0	1.03
26	58.7	93.2	0.95
27	65.5	103.1	0.85
28	22.7	122.5	1.87
29	24.0	130.6	1.68
30	24.9	137.1	1.62

Results may vary depending on materials and methods used.

**1.800.267.2796**

**info@precisionbiologic.com**